



## IMD-A<sup>®</sup> Series

Rapid Biologic Detection System for Air

**The IMD-A family of Real-Time, Rapid Biologic Detection Systems for Environmental Air Monitoring can continuously guard aseptic pharmaceutical environments to instantly warn when microbes are present**

- *Real-time operation*
- *Measures particle size and biologic status simultaneously*
- *Immediate results without staining or reagents*
- *Continuous monitoring, or by specified sample volume*
- *Visual status indicator*
- *Built-in system self test*
- *USP<1223> compliant*
- *21 CFR Part 11 compliant*
- *Stainless steel enclosure*
- *Sanitary flange for easy removal and re-installation*
- *Flexible communications interface for networking*



[www.metrovalidpharma.com](http://www.metrovalidpharma.com)


Résidence Gimmo, Immeuble N°03 Bab Ezzouar, Alger - 16042 - Algérie  
Tél. : +213 23 88 46 43 à 45 /Fax : +213 23 88 46 42

*All information contained in this data sheet is being presented for informational purposes only, and is subject to change without notice.*

## BioVigilant's IMD-A rapid biologic detection system revolutionizes environmental air monitoring for microbial contamination.

Using patented optical technology, the IMD-A system is capable of simultaneously detecting both size and biologic fluorescence of airborne particulates. The technology offers an eloquent yet robust method for determining the microbial content of an aerosol sample on a particle-by-particle basis, thereby providing a continuous, real-time analysis of environmental bioburden.

The IMD-A system requires no reagents and minimal human operator involvement, so it is particularly well-suited for rigorous environmental monitoring applications where time, cost, and cleanliness are of greatest concern.




Microbes contain metabolites or chemical compounds such as NADH, riboflavin, and dipicolinic acid (DPA). When excited by a certain wavelength of light, these compounds emit an intrinsic fluorescence signal that serves as a marker for biological content of the sample.

### *Bacillus Atrophaeus under Confocal Microscope*

When a particle enters the IMD-A system's interrogation zone, it intersects the beam of a single laser light source.

The scattering of the laser light determines the size of each individual particle. At the same time, the presence or absence of intrinsic fluorescence determines whether the particle is biologic or inert.




### *IMD-A's Optical Subsystem*

[www.metrovalidpharma.com](http://www.metrovalidpharma.com)

*The supplied sanitary flange clamp at the inlet of the IMD-A system allows for easy connection to tubing, and quick disconnection during cleaning.*



*The IMD-A system also features a video option that allows simultaneous video capture of the environment being measured for precise correlation of physical events with detected microbial or particulate excursions.*



## Communicating with IMD-A System

Cables supplied with these peripherals

Option 1 (Recommended): Connect IMD-A locally to standard monitor and keyboard via IMD-A USB connections.

Crossover cable

Option 2: Connect IMD-A locally to external computer via crossover cable. (Requires use of a remote desktop application to access PharmaMaster software application.)

Ethernet cable

Option 3: Connect IMD-A to remote computer via company network via Ethernet. (Requires use of a remote desktop application to access PharmaMaster software application.)

## The IMD-A system is supplied with PharmaMaster software.

### Instant-Read Interface

The system's sophisticated PharmaMaster software offers different levels of functionality based on user need and authority. For environmental control, the software features graphic displays that visually communicate custom-defined limits of particle contamination. A green light indicates that operation is safely within defined limits, yellow warns that an Alert level has been exceeded, and red indicates that an Action level has been exceeded. Alert and Action levels can trigger alarms, if desired.

### User-Defined Alert Protocols

PharmaMaster software enables users to create a library of alert protocols with sensitivity levels to match different applications and environments.

### Results Displayed by Particle Size


With PharmaMaster software (as shown in the display at right), users can view independently controlled line graphs to display results for all particles greater than 0.5 µm, all particles larger than 5.0 µm, and all biologic particles. An alarm alerts managers when a user-defined level is exceeded.

### Five-Channel Sorting

With the PharmaMaster software, users can choose to review biologic particles in five data bins, from 0.5 µm to >10 µm.

### Two Operational Modes for Application Flexibility

The PharmaMaster software offers both a continuous detection mode and a point sample mode for a sampling a specified volume or time duration, allowing users to determine the method of environmental sampling that best fits their application.



Graphical User Interface

## Ordering Information

### Systems

- IMD-A 300: IMD-A with 1.15 LPM (±5%) sample flow rate
- IMD-A 350: IMD-A with 28.3 LPM (±5%) sample flow rate utilizing virtual impaction concentration method

### Accessories

- A01000 Panasonic Toughbook computer
- A02000 Monitor, DVI-compatible
- A02010 Sealed keyboard with mouse
- A02020 Crossover cable, 6ft.
- A03000 Sanitary flange kit
- A03001 Isokinetic probe for IMD-A 300
- A03002 Isokinetic probe with stand for IMD-A 300
- A03003 Isokinetic probe for IMD-A 350
- A03004 Isokinetic probe with stand for IMD-A 350
- A04001 Bev-A-Line XX tubing, 3/8", 100ft. (recommended)
- A04002 Bev-A-Line XX tubing, 1/2", 100ft.
- A04003 Disposable HEPA filter with Polypropylene housing

## Services and Support

### System Verification Templates

Pre-printed documentation sets for system functional design specification (DQ), installation qualification (IQ), and operational qualification (OQ) are available for completion without the need to recreate these documents.

### Onsite Support Services

BioVigilant's highly trained applications and service personnel can streamline maintenance, implementation and validation work. Our experts can assist with any of these activities:

- Calibration
- Preventive maintenance
- Equipment installation
- Test plan development
- Data analysis
- Operator training
- Installation Qualification (IQ)
- Operational Qualification (OQ)
- Performance Qualification (PQ)

Contact your sales person to learn more.

## Technical Specifications

Detection Method:	Optical, Mie Scatter (particle size) Optical, Auto-fluorescence (biologic status)	Data storage:	Peripheral hard drive (≥80GB)
Biologic Detection Capability:	Vegetative cells, bacterial spores, fungi	Communications output options:	(1) Ethernet 10/100-base T (4) USB 2.0 ports (1) Digital Video Interface (DVI)
Biologic Counting Efficiency:	Meets USP<1223> and EP 5.1.6	User interface:	Compatible with standard DVI/SVGA monitor, mouse, and keyboard.  Remote access and control via network
Incubation time:	None required	Software:	Windows® XP Operating System, PharmaMaster® Bio-Detection Software, VNC® Server Enterprise Edition, and SQL Server® Express Edition
Particle Size Range:	0.5µm to ≥10µm	Data export:	ASCII, CSV
Particle Sizing Bins:	Grouped by size into 5 bins: 0.5µm to 1µm 1µm to 5µm 5µm to 7µm 7µm to 10µm ≥10 µm	Languages supported:	English, Japanese
Sampling Mode:	Continuous or Point sampling for specified duration or volume	Regulatory compliance:	Electromagnetic compatibility 2004/108/EEC CE/GS Mark RoHS 2002/95/EC WEEE 2002/96/EC 21 CFR Part 11 FDA Type V Drug Master File
Particle Threshold:	0.5µm	Standards met:	USP<1223> Validation of Alternative Microbiological Methods EP 5.1.6 Alternative Methods for Control of Microbiological Quality Manufactured in accordance with ISO 13485
System Dimensions:	9" W x 18 1/8" H x 24 1/4" D	Accessories included:	Bev-A-Line XX tubing (3/8" ID), (8) feet Sanitary flange kit for 3/8" tubing Inlet HEPA filter Web camera Inlet and outlet caps User manual
Weight:	54.0 lbs (24.6kg) (IMD-A 300); 55.5 lbs (25.2 kg) (IMD-A 350)		
Enclosure material:	304 stainless steel, Nickel-plated aluminum		
External chemical resistance:	Isopropyl alcohol (91%), chlorinated solution (2500ppm), hydrogen peroxide (3%)		
Operating temperature range:	15-30 degrees Celsius		
Operating Humidity Range:	10-80 % RH Non-condensing @ 15-25 degrees Celsius, 10-60% RH Non-condensing @ 25-30 degrees Celsius		
Altitude:	< 6000 feet (1830m)		
Power Requirements:	100-240 VAC, 50-60 Hz		
Power consumption:	90 Watts, max.		
Air Outlet filtering:	HEPA filtering (Internal)		
Light source:	405nm laser		
Detectors:	Photodiode, Photo Multiplier Tube		
CPU:	Embedded Intel® dual core 1.5GHz, 4MB L2 cache		

BioVigilant, the BioVigilant Logo, IMD, IMD-A, PharmaMaster, and the term Instantaneous Microbial Detection are trademarks or registered trademarks of Azbil BioVigilant, Inc., in the United States and/or other countries. Other names or brands may be the property of others.  
LI005 December 4, 2013